



Milk



Unlock next-level testing with Bioeasy range

Comprehensive, accurate, reliable and simple testing for antibiotic, aflatoxin M1, and adulterant residues in raw milk — powered by Novonesis

BioeasyTM

Test Brochure

Test kits /Antibiotics group	Incubation time	Incubation temperature	Type	Test protocol	Antibiotics
Beta-lactams test - YRM1001-2C	2 + 3 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams
2IN1 BT (EU) - YRM1007-40	3 + 5 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines
2IN1 BTCef (EU) - YRM1008-40	3 + 6 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines Cephalexin Chloramphenicol
4IN1 TCCB - YRM1126	1 + 5 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines Cephalexin Ceftiofur
3IN1 BST - YRM1026	3 + 6 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines Sulfonamides
4IN1 BSCT - YRM1039C	5 + 7 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines Streptomycin Chloramphenicol
4IN1 BCTC - YRM1039CQ	3 + 7 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines Chloramphenicol
4IN1 BTSQ - YRM1059	3 + 6 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines Sulfonamides Quinolones
2IN1 BT (EU) - YRM1007-40T	3 + 3 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines
4IN1 BSCT - YRM1039C-10	3 + 7 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines Streptomycin Chloramphenicol
2IN1 BT (RU) - YRM1008-1C	3 + 5 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines
Beta-lactams & Tetracyclines - YRM1096	3 + 4 min	40°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines

The information provided herein is for general information purposes only. All information is provided in good faith. There is no warranty being extended and no liability whatsoever (including without limitation, any direct or indirect damages for lost profits or business interruption) is accepted as to its accuracy, completeness, correctness, non-infringement, merchantability, or fitness for a particular purpose. The product(s) may be covered by pending or issued patents, registered or unregistered trademarks, or similar intellectual property rights. Copyright © Novonesis Group, by ownership or license. All rights reserved.

BioeasyTM

Test Brochure

Test kits /Antibiotics group	Incubation time	Incubation temperature	Type	Test protocol	Antibiotics
2IN1 BT (EU) - AEMRT1007-40 (ambient)	3 + 5 min	18-25°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines
2IN1 BTCef(EU) - AEMRT1008-40 (ambient)	3 + 6 min	18-25°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines
3IN1 BST - AEMRT1026 (ambient)	3 + 6 min	18-25°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines Sulfonamides
4IN1 BTSQ - YRM1059RT (ambient)	3 + 6 min	18-25°C	A lateral-flow standard strip test	2 step	Beta-lactams Tetracyclines Sulfonamides

Test kits /Antibiotics group	Incubation time	Incubation temperature	Type	Test protocol	Antibiotics
Florfenicol & Thiamphenicol - YRM1054-1	3 + 3 min	40°C	A lateral-flow standard strip test	2 step	Phenicol
AMINO 3IN1 - YRM1034	3 + 6 min	40°C	A lateral-flow standard strip test	2 step	Aminoglycosides
Gentamicin - YRM1007-100	3 + 5 min	40°C	A lateral-flow standard strip test	2 step	Gentamicin
MACRO 3IN1 - YRM1022	3 + 7 min	40°C	A lateral-flow standard strip test	2 step	Macrolides
4IN1 CGSN - YRM1098	3 + 6 min	40°C	A lateral-flow standard strip test	2 step	Chloramphenicol Aminoglycosides
Fluroquinolones - YRM1006-20	3 + 3 min	40°C	A lateral-flow standard strip test	2 step	Quinolones
Aminoglycosides - YRM1080	3 + 7 min	40°C	A lateral-flow standard strip test	2 step	Aminoglycosides
Chloramphenicol - YRM1010-8	3 + 5 min	40°C	A lateral-flow standard strip test	2 step	Chloramphenicol
2IN1 QS - YRM1016A	3 + 5 min	40°C	A lateral-flow standard strip test	2 step	Sulfonamides Quinolones

The information provided herein is for general information purposes only. All information is provided in good faith. There is no warranty being extended and no liability whatsoever (including without limitation, any direct or indirect damages for lost profits or business interruption) is accepted as to its accuracy, completeness, correctness, non-infringement, merchantability, or fitness for a particular purpose. The product(s) may be covered by pending or issued patents, registered or unregistered trademarks, or similar intellectual property rights. Copyright © Novonesis Group, by ownership or license. All rights reserved.

BioeasyTM

Test Brochure

Test kits /Antibiotics group	Incubation time	Incubation temperature	Type	Test protocol	Antibiotics
Aflatoxin M1 (15-150 ppt) - YRM1004D-2	2 + 8 min	40°C	A lateral-flow standard strip test	2 step	Aflatoxin M1
Aflatoxin M1 (300-400 ppt) - YRM1004-3	3 + 4 min	40°C	A lateral-flow standard strip test	2 step	Aflatoxin M1
Aflatoxin M1 (40-50 ppt) - YRM1004-0.05B	3 + 6 min	40°C	A lateral-flow standard strip test	2 step	Aflatoxin M1
Aflatoxin M1 (100-1000 ppt) - YRM1004D-4	2 + 8 min	40°C	A lateral-flow standard strip test	2 step	Aflatoxin M1
Aflatoxin M1 (150-600 ppt) - YRM1004D-1	2 + 8 min	40°C	A lateral-flow standard strip test		Aflatoxin M1

Test kits /Antibiotics group	Incubation time	Incubation temperature	Type	Test protocol	Antibiotics
					Sulfonamides Fluoroquinolones Aminoglycosides Phenicol Macrolides
16IN1 / YRM1112	3 + 7 min	40°C	A lateral-flow standard strip test	2 step	
Melamine (80-100) - YRM1003-100	3 + 3 min	40°C	A lateral-flow standard strip test	2 step	Melamine
Melamine (30-50) - YRM1003-50	3 + 3 min	40°C	A lateral-flow standard strip test	2 step	Melamine
Whey adulteration in milk - YREC1003	3 + 3 min	40°C	A lateral-flow standard strip test	2 step	Cow whey
Whey adulteration in milk - YREC1003-3	3 + 3 min	40°C	A lateral-flow standard strip test	2 step	Cow whey
Cow milk adulteration in goat & sheep milk - YREC1001-4-Y	3 + 3 min	40°C	A lateral-flow standard strip test	2 step	Cow milk
Cow milk adulteration in goat & sheep milk - YREC1001-5-Y	3 + 3 min	40°C	A lateral-flow standard strip test	2 step	Cow milk

The information provided herein is for general information purposes only. All information is provided in good faith. There is no warranty being extended and no liability whatsoever (including without limitation, any direct or indirect damages for lost profits or business interruption) is accepted as to its accuracy, completeness, correctness, non-infringement, merchantability, or fitness for a particular purpose. The product(s) may be covered by pending or issued patents, registered or unregistered trademarks, or similar intellectual property rights. Copyright © Novonesis Group, by ownership or license. All rights reserved.